[®]*Dentag*

00

00

2

.A0. Manico Anatomico

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Anatomic Handle

Obtained from a stainless - steel rod with a diameter of mm 5.5 or 6 into which the tips of the instrument has been inserted.

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Serrated Hollow Round Handle mm 8

Obtained from a stainless - steel tube with a diameter of mm 8 into which the tips of the instrument have been inserted. This type of handle is particularly light.

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Round mm 10 Serrated Hollow Handle

Obtained from a stainless - steel tube with a diameter of mm 10 into which the tips of the instrument have been inserted. This type of handle is particularly light and comfortable.

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One piece instrument

Instrument obtained from a single, round, stainless - steel rod.

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DenTag Handle mm 6

Smooth handle, obtained from a stainless – steel rod with a diameter of mm 6, into which the tip of the instrument has been inserted.

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DenTag Handle mm 8

Smooth handle, obtained from a stainless – steel rod with a diameter of mm 8, into which the tip of the instrument has been inserted. This handle is particularly light.

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DenTag Handle mm 10

Smooth handle, obtained from a stainless – steel rod with a diameter of mm 10, into which the tip of the instrument has been inserted. This type of handle is particularly light and comfortable.

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Octagonal handle

Instrument obtained from a single, octagonal, drawn piece of stainless – steel with a diameter of mm 5.5.

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Round Aluminium Handle mm 8

Obtained from an aluminium rod with a diameter of mm 8 into which the tips of the instrument have been inserted. This type of handle is available in various colours to enable the instrument to be identified at first glance. It is also extremely light.

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03/04 COL

The **DenTag** code is made up of three distinct parts that are separated by a period.

1) The first part, which consists of four digits, identifies the type of instrument (scaler, probe, elevator, etc.).

3400

Gracey Curettes

2) The second part, which consists of letters and digits, identifies the type of handle used on that particular model of instrument.

DT08

DenTag hollow round handle mm 8 fig. **3/4 red**

3) The third part, which generally consists of digits or letters, is the international identification code for the instrument. This code describes its length and colour, also contains other identifying information.

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Warnings

Each and every precision instrument manufactured by DenTag for use in dental surgeries is produced and monitored with the utmost care by highly-skilled technicians. In order to guarantee the maximum efficiency and service life of these instruments, it is necessary to use them with care and solely for the purpose for which they are designed.

Special care must be taken when cleaning and sterilizing the instruments:

- only specific products should be used; any chemical-based products must be used in the recommended concentrations and according to the manufacturer's instructions;
- temperatures must not exceed 175 °C - 347°F.

It is important to underline that certain types of damage (corrosion, rust spots and similar) may be caused to the instrument by simple contact. Failure to observe these fundamental rules may result in irreparable damage to the instrument and would render null and void any warranty and/or claim against DenTag.

Important: The instruments are supplied unsterilized, and must be inspected, cleaned and sterilized before use. After using the instrument, it must be inspected for corrosion, rust spots, broken, loose or bent parts, wear of the blade edge or other damage. Any instrument which is not in perfect condition must immediately be removed from service and conserved separately. Disposable instruments must never under any circumstances be sterilized and re-utilized, but must be removed from service and separated immediately after use.

Cleaning and maintenance

Brief recommendations concerning the cleaning and sterilization of our instruments, which are supplied UNSTERILIZED, and must always be washed and sterilized before use.

The instruments used must be rinsed under hot running water immediately after use, and any organic residue must be carefully removed. The instruments must be cleaned immediately after rinsing. If this is not possible, they must be immersed in a solution of water and a neutral pH detergent.

The instruments may be cleaned using a washing machine, by immersion in an ultrasonic tank, or by hand. The following general guidelines should be followed:

Cleaning in a washing machine

Follow the manufacturer's instructions. Ensure that hinged and jointed instruments are correctly lubricated before the sterilization cycle.

Cleaning in an ultrasonic tank

Immerse the instruments completely in the washing solution and observe the times recommended by the manufacturer (generally 5 – 10 minutes).

Ensure that instruments featuring hinges or joints are placed in the open position in the tank. Instruments with blades – scissors, scalpels, curettes etc. – must never be placed in contact with each other or other instruments. Instruments made of different metals should not be washed together. Do not use the washing solution for too many washing cycles. Follow the manufacturer's instructions and replace the solution frequently. At the end of the washing cycle, eliminate any washing solution residue by rinsing thoroughly under running water.

Manual cleaning

If it is not possible to clean the instruments in a washing machine or ultrasonic tank – which are the strongly recommended methods – the instruments may be washed by hand. In this case:

- use brushes in rigid nylon, plastic or similar;
- never use metal brushes or steel wool;
- always use neutral pH detergents.

When the instruments are perfectly clean and rinsed, they must be carefully dried. All hinged or jointed instruments must be lubricated using lubricants designed specifically for surgical instruments. Never use household or industrial oils and lubricants. Check the condition of each instrument and that it functions correctly.

Only now can the instruments proceed to the sterilization phase, which may be carried out either in batches or individually.

High-temperature sterilization in batches

The instruments are placed in the special trays and/or baskets and placed in the autoclave. All hinged or jointed instruments must be introduced in the open position. (A closed instrument will never be sterile, since the steam cannot penetrate between the parts in contact).

High-temperature sterilization of individual instruments

The instruments are placed individually in special disposable bags. Hinged or jointed instruments must be introduced into the bags in the open position.

Cold sterilization

The instruments are sterilized by immersion in a sterilization solution for a specific period of time as specified by the manufacturer of the solution. The manufacturer's instructions must be followed carefully.

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Guarantee declaration

- All the products are bearing their own item code
- Every item of our instruments is manufactured using the same material
- All our instruments are bearing our trade mark, the CE mark, their particular code and the number of lotto for tracking use. (as the CE rules want)
- Our instruments can be sterilized in sterilizing machine, in autoclave and with chemical solutions (not sodium hypochlorite or similar)
- The raw materials we use respect the UNI EN ISO 7153-1:2016
- The range of hardness our stainless steel is between 40-52 Rockwell- (According to DIN 9001- DIN 58298 Rockwell Hardness Standards)
- The quantities of impurities included on the stainless steel are as from the Standard Characteristics.
- Our instruments are manufactured according to UNI CEI EN ISO 13485:2012 Quality Systems – Medical devices, and CE rules
- Our instruments can be sterilized in 134 C° autoclaves which works with steam and also with gas
- Our instruments are covered by guarantee against corrosion, oxidation and fractures caused by defects of manufacture for 5 (five) years. This guarantee do not cover logically every improper use
- Test of corrosion made by us; Boiling in distilled water for 1 hour and drying in air. This 2 or 3 times. After we control with optical instrument. The only positive result after this test is the absolute absence of rust.